

## Problem Solving

Mechanical Problem	Possible Cause	Action
<b>Bag will NOT infuse</b>	Clamp closed on transfer set Blue clamp on new bag line Kink in tubing Cones not completely broken Internal organs blocking flow	Open clamp Remove clamp Straighten tubing Rebreak cones Change your position, cough, raise bag
<b>No Drain</b>	Clamp closed on transfer set Blue clamp on drain closed Cones not completely broken Kink in tubing Fibrin in tubing Constipation Internal organs blocking flow	Open Clamp Open Clamp Rebreak Cones Straighten tubing "Milk" tubing Follow MD's orders Change your position
<b>Contamination</b>	Touch y set end Transfer set dropped or touched closed	Throw bag away, start over Make sure clamp is closed <b>DO NOT</b> do exchange Apply new disconnect clamp, call clinic, and come for transfer set change